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MORVABON.COM

Calci Mor

Light-cured Resin Base
Cavity Liner

manufacturer

Pardis Morvabon Pharmed Co.:

Head Office: Unit 154, 7nd floor, c block,
Borjsaz Bd, Azadi Street, Tehran, Province, Iran
Tel: 0098 21 6658 1438
Fax: 0098 21 6690 8819

- ISO13485 : 2016 , US cert, IAS
- ISO 9001: 2015, TÜV certification
- ISO 14001: 2015, TÜV certification
- ISO 18001: 2007, TÜV certification
- ISO 10002: 2014, TÜV certification
- ISO 10004: 2012, TÜV certification

MEDICAL IPDA
EQUIPMENT
DEPARTMENT

PRODUCTION
LICENSE NO.
4 9 8 2 5 8 7 2
IMD: 69400836

ISO 13485: 2016

Calci Mor

Light-cured Resin Base Cavity Liner

Calci Mor is a resin based, light-cured, radiopaque cavity liner, designed for use in direct and indirect pulp capping and as a protective base/liner under dental adhesives, varnishes, filling materials, cements, and other materials.

Composition:

resins and Barium Sulfate.

Uses:

- ▶ Protect pulp of the tooth
- ▶ Lining under all filling materials
- ▶ Protection when applying the total etch technique

Benefits:

- ▶ Superior Physical Properties
- ▶ Low Water Solubility
- ▶ Radiopacity
- ▶ Light curing type cavity liner cure extremely hard
- ▶ Easy dispensing

Instruction for use:

- ▶ Prepare the cavity in a conventional manner.
- ▶ Rinse the cavity and dry
- ▶ Apply directly to the cavity with the syringe. Place a new lure lock.
- ▶ Remove any material inadvertently applied over enamel margins.
- ▶ Cure 30 seconds, keeping the tip as close to the cured surface as possible. (with halogen or LED light with a minimum intensity of 400 mW/cm² in the 400-500 nm range.)
- ▶ Apply a second layer if indicated.
- ▶ Cure for 30 seconds.

Storage:

- ▶ **Calci Mor** material containers should be stored with the original cap only.
- ▶ kept out of direct sunlight and it would be better to store in a well-ventilated place at refrigerated temperature between 2°C and 8°C when not in use.

- ▶ Allow material to reach room temperature prior to use.
- ▶ Do not freeze.
- ▶ Protect from moisture.
- ▶ Do not use after expiration date.
- ▶ Keep out of reach of children.

Precautions:

- ▶ Only for dental use.
- ▶ In case of consumption drink copious amount of water. Do not induce vomiting. Seek medical advice.
- ▶ Wear suitable protective clothing, eyewear, and gloves.
- ▶ Is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins.
- ▶ Do not place on enamel, cavity margins or leave exposed to oral environment.
- ▶ The success of pulp capping procedures may be influenced by preparation contamination. Use of a rubber dam or adequate isolation is recommended.
- ▶ If a thickness of material greater than 1mm is desired, it should be applied and cured incrementally. As the material increases in thickness, the shade becomes more intense (white).
- ▶ Tooth size, shape, dentinal radiopacity, exposure conditions, thickness of the Calci-Mor and overlying restorative material may affect the ability to radiographically detect the material.
- ▶ The container should be tightly closed immediately after use.
- ▶ Use only in well ventilated areas.

					
Consult instructions for use	Indicates the batch code	Catalogue number	Indicates the date after which the medical device not to be used.	Indicates a medical device that needs protection from light sources.	Indicates a medical device that needs to be protected from moisture.
					
Do not use if Package is Damaged	Indicates the medical device manufacturer (Name & Address)	Indicates a medical device that has not been subjected to a sterilization process.	Indicates the device intended for use by professional person only.	Indicates that the medical device can cause discomfort and reddening, irritation, or blistering of the skin	Indicates the Temperature Limit

Produced by Morvabon Co.
 www.morvabon.com
info@morvabon.com
 Made in IRAN