



Instagram



Moheito



YouTube



Youtube-EDA Channel



Facebook





Changzhou Sifary Medical Technology Co., Ltd

Add: NO.99, Qingyang Road, Xuejia County, Xinbei District, Changzhou City, 213000 Jiangsu, China Email: ivy@sifary.com Web: www.eighteeth.com





# CATALOG

HyperLight Portable X-ray unit	01-04
NanoPix Intraoral X-Ray Sensor	05-08
Helios 600 Intraoral Scanner	09-14
AirPex Wireless Apex Locator	15-16
E-CONNECT & E-PEX Cordless Endomotor & Apex Locator	17-22
E-connect S Endomotor With Built-in Apex Locator	23-24
E-value Brushless Endo Motor	25-26
E-xtreme Make Essence to Extreme	27-28
Fast-Pack & Fast-Fill Cordless Obturation System	29-36
Ultra X Cordless Ultrasonic Activator	37-38
UltraMint & UltraMint Pro All-round Ultrasonic Scaler	39-44
CuringPen 4 LEDs Curing Light	45-46







## **Engineered for accurate images**

- Clinical Imagery








- Sharp, high resolution images with a *0.4mm* focal spot

Source focal spot Source focal spot





Effect of focal spot size



## Designed with safety in mind



Built-in leaded double shield to protect the operator from leakage radiation



Locking system to prevent unintentional exposure

## **Specifications**

Model	HyperLight	Weight	1.8 kg
Focal Spot	0.4 mm	Battery	2500 mAh x 4
Tube Voltage	65 kV	Exposure Time Range	0.02-2 Sec
Tube Current	2.5 mA	Line voltage	100V-240V, 50Hz/60Hz


HyperLight 02







## Simplify your digital workflow


- Compact size and lightweight can make you operate it by one hand only
- Intuitive user interface












04









1/10 Cost per patient compared to film 3s To complete entire imaging process 50000+ Exposures with no quality loss



## Leading technologies: The 1st time in dentistry



#### AED Technology (Automatic Exposure Detection)

- from chest radiography
- Compatible with almost all DC and AC X-Ray generators
- Digital images comes automatically when your X-Ray is switched on



#### **APS CMOS Technology**

- from aerospace industry
- High signal-to-noise ratio with maximum details
- Long life and low power consumption



#### Direct Deposition CSI Technology

- from chest radiography
- Low dose radiation for safety
- Less exposure time but higher resolution

All-in-one
The most economical X-Ray sensor in the world



#### Unique design



#### Ergonomic design

Cupped cable connection ensures comfortable finger positioning



#### Robust cable connection

70000+ Bendings RA test passed



#### Ultra-thin body design

4.4mm - the thinnest X-Ray sensor

NanoPix

06

## Easy-of-Use software with only

## 4 steps

Step 1. Choose 'Patient'



Step 3. View Image



Step 2. Image 'Acquisition'



Step 4. Generate 'Report'





## The Complete Image Solution:

excellent combination of DR Pad & its base




pre-installed software



multiple USB ports for diverse needs



build an instant connect with users



## Specifications(available in 2 sizes)

	NanoPix 1	NanoPix 2
Detector Technology	APS CMOS	APS CMOS
Scintillator	Csl	Csl
Pixel Matrix	1000*1500	1300*1800
Theoretical Resolution	25lp/mm	25lp/mm
True Resolution	20lp/mm	20lp/mm
Active Area (mm)	20*30	26*36
Mechanical Size(mm)	25.4*36.8*4.4	30.4*41.9*4.4
AD Conversion(bit)	16	16
Data Interface	USB 2.0	USB 2.0
Protection	IP68	IP68





# Start your smooth digital journey with Helios 600!

#### Why bother?

High accuracy with automatic and adaptive calibration, manual calibration is not necessary any more.

#### Ready for your digital experience-

One of the most accurate and fastest scanners on the market.

#### Ready for your hand-

Small enough to operate it easily, though it achieves a 16mm x 14mm high-definition field of view.





#### Lighte

198 grams, one of the most lightweight scanners on the market, well-balanced, ergonomically designed, and easy to carry.



#### Faster

With Helios 600 you can take a full arch scan in less than 30 seconds.




#### Smarter

Al technology filters the soft and moving tissues, offering users a smooth scanning experience and generating a natural-looking treatment outcomes.




#### More accurate

As one of the most accurate scanners on the market, Helios 600 has reached the state-of-the-art accuracy level- the average deviation of a in-vivo full arch scan is around 20µm.







## Plug & Scan



Smart Reminder



Air mouse



Automatic software update



Timely after-sales service



Cloud Storage



Voice Memo



Online 3D viewer

## **Specification**

#### Helios 600

Accuracy	Full-arch 20.0μm		
C	16 x 14mm (L)		
Scan area	12 x 12mm (S)		
	Dimension	219 x 46 x 36mm	
Handpiece	Weight	198g	
	Capture button	Two sided	
Tim	Tip size (front part)	22 x 18mm (L)	
Tip		18.6 x 16.7mm (S)	
	Length	1.8m	
Connectivity	Туре	USB 3.0	
C-III	Accuracy calibration	Calibration-free	
Calibration	Color calibration	Calibration-free	
Special mode	Remote control mode		

#### Recommended System Requirements

Processor	Intel®CoreTM i7 9th Generation, base frequency 2.6 GHz (or better)
Memory	16GB (or more) DDR4, frequency 2666 MHz (or better)
Hard Drive	512GB (or more) SSD
Graphics card	NVIDIA®GeForce®GTX 1650 (or better)
Display	15.6" FHD (1920 x 1080) (or more)
Others	USB 3.0 port
Operating system	Windows 10 Pro
Optional	Touch screen

Helios 600




## AirPex

## Wireless Apex Locator

International Patent No.:PCT/CN2019/110429







## Feel the Charm of Wireless Charging

Lay it down , Charge it up Charging is as simple as placing the AirPex on the charging base



## **True Wireless Apex Locator**

Get rid of the wire to the fullest Use without magnet wire Accuracy through wirelessness



# 88 E-CONNECT & E-PEX Cordless Endomotor & Apex Locator

E-CONNECT & E-PEX





**Torque Calibration** 



Powerful 1500 mAh Lithium Battery



Digital OLED Display



Integration Option



Low Noise Motor



Screen Cast / Mirroring



Ambidextrous \* Display

\*User can rotate the Endomotor Screen towards left or right \*In connect function Endomotor mirrors Apex Locator's screen



#### Five Reciprocation Modes (MO)



Reverse Reciprocation-One Mode Counter clockwise & clockwise movement

Fwd 30° - Rev 150° For reciprocating files



Fwd 150° - Rev 30° Fwd 180° - Rev 30° Fwd 210° - Rev 30° Fwd 250° - Rev 30° For all continuous rotation files



#### Nine Editable Memory Programs (M1 to M9)

















#### One Piece Contra Angle Head



340° Adjustable Handpiece



Detachable Head



Miniature Head



Autoclavable

19 20 E-CONNECT E-CONNECT







## Integrate Endomotor with Apex Locator (Connect Function)









#### **Smart Accurate and Reliable Apex Locator**

- Advanced Multi frequency Technology
- High Precision in Wet & Dry Canals
- Large 3.5" LCD Display
- Automatic Calibration
- Compact Size & Stable Design

#### Screen Display & Its Meaning

	Interface	Interface	Interface	Interface
Display	3 - (v)	3 2 4 (v)	3 2 4 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	3 2 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
Value	12	09	00	-3
Meaning	The corresponding scale value is 2, indicating the distance from the apex is about 2mm	It indicate that the distance from the apex is about 1.5mm	It indicates that the endo file reaches the apex	It indicates that the root canal file has exceeded the apex
Alarm	Didi alarm sound for a long time interval	Didi alarm sound interval becomes shorter	Uninterrupted long beep	Very rapid Didi alarm sound

21 E-PEX 22



## E-connect S

Endomotor With Built-in Apex Locator





Torque Calibration
For precise torque output



OLED Screen
Digital OLED Display



Ambidextrous Display
Screen can be rotated as per
Left or right hand use




Powerful Battery Heavy duty 1500mAh Lithium battery





Apex Locator Bar



#### APEX LOCATOR

- Apical control with rotary instruments
   Position of file is displayed on screen and the motor rotates in reverse direction when file reaches pre-set apical working length, to prevent over instrumentation
- Manual canal measurement
   Manual file clip can be used to measure canals precisely with hand files
- Multi frequency technology with high accuracy & precision
- Apex locator flash bar with zoom-in function

## **Working Modes:**

**FWD** - Clockwise Continuous Rotation

**REV** - Anti-clockwise Continuous Rotation

**REC** - Reciprocation Movement

First-in-class cordless open motor with adjustable angles of rotation ranging from 30° to 370° for reciprocation as per user's need  $\,$ 

First-in-class feature of speed adjustment in reciprocation mode

#### ATC - Adaptive Torque Control

ATC manages the file stress in a smart way. When torque limit is reached, instead of going into reverse mode, motor reciprocates till the stress is reduced. Once back to below pre-set torque level, file will again rotate in continuous rotation

#### Integrated mode features:

- Auto Start & Stop Motor starts working when the file enters inside root canal, and stops when withdrawn
- Apical Reverse When the file reaches pre-set apical reverse point, user can define whether the file starts moving in reverse direction
- Apical Slow Down To preserve the cemento-dentinal junction and avoid file separation, speed of motor slows down when file is in apical region

#### **Specifications:**

Rotational Speed Range : 120-1000 rpm Torque Range : 0.5-4.0 N.cm Memory Programs : 11 (M0 to M10) Battery Capacity : 1500 mAh

Eighteeth




## E-value

#### **Brushless Endo Motor**

Eighteeth

value



25

Speed Range : 120-1000 rpm

Torque Range : R.L-4.0 N.cm

Memory Programs : 10 (M1 to M10)

Battery Capacity : 1900 mAh



OLED Color Screen

More comfortable visual experience



Built-in File System Switch file systems quickly



360° Rotatable 360° Rotatable Contra Angle



Brushless Motor Reliable, Durable, Efficient



Long Battery Life
Heavy capacity
1900mAh lithium battery




Mini Contra Angle Advanced 1:1 contra-angle One piece design





#### Built-in File System




Built-in File System

Long Press S to enter the Mode Setting interface

E-FLEX BLUE M1 300rpm
Fwd 3.0N · cm

Choose the Endo File being used and press S

Note: the doctor could adjust the speed, torque limit and reciprocating angle according to the needs of the case.



#### Working Modes:

• Fwd - Forward Rotation

• Rev - Reverse Rotation

REC - Reciprocation Mode

Doctor could set the reciprocating angle and speed level according to actual needs

- Speed Level range from Level 1 to Level 5
- -Angles range from 30° to 370°
- ATC Adaptive Torque Control Mode

After reaching the preset torque limit, motor will reciprocate according to the set angle to improve the efficiency of root canal enlargement







26

E-value



## E-xtreme

Make Essence to Extreme

Smaller | Lighter | Quieter than ever





360° contra-angle rotation



Slim and compact design for maximum comfort and mobility



Miniature contra-angle for better visibility & easier pos t erior access

Miniature contra-angle Super Mini head & thin neck



- Normal head & neck





Eighteeth innovative motor, cutting down the noise from inside, allowing you to experience the tranquility.



Less than 100 grams, reducing hand fatigue during treatment, providing the ultimate E-xtreme experience.





Reciprocation function





Built-in file system

## **Specifications**

Model : E-xtreme : 99 g Weight(grams)

Power supply

: 100 V-240 V

Torque range : 0.5 N.cm-4.0 N.cm Output

Frequency

: 50 Hz / 60 Hz : DC 5 V/1A

27

E-xtreme

E-xtreme

28





## Fast-Fill

Obturation Fill Handpiece

**P**Eighteeth

Fast-Fill



#### Fast-Fill Tips

- Autoclavable silver tips for multiple usages
- Silver tips maintain the thermoplastic temperature of GP and offers smooth flow inside canal
- 360 ° rotating tip for easy filling at all angles
- Extra long needle offers convenience and visibility
- Easy to maintain and replace

#### Motorized GP Flow



- Provides consistent flow of GP at desirable rate
- Precludes voids and makes obturation homogeneous
- Eliminates hand fatigue that resulted from earlier gun type devices

#### Comprehensive Display



- Three customizable memory programs T1 to T3
- Temperature range 100 °C 200 °C
- GP flow rate Low, Medium, High
- Reserve GP level indicator

## Intuitive Design - Sleek and ergon

- Sleek and ergonomic hand-piece
- Simple push button to dispense GP
- LED indicator reflects when the device is ready for use

#### Powerful Battery



- Highly efficient 2600mAh Lithium battery
- Direct Plug & Play In case of low battery, use the device directly with adaptor

#### Specifications

Heating Time	20 seconds to 200 °C
Charging Time	4 hours to complete charge
Hexagonal Wrench	Use the wrench to bend the tip, rotate the tip and clean the device



23 Gauge Large



25 Gauge Small





#### **Temperature Range**



- Wide temperature range of 90° to 250°
- Five memory programs to customize temperature and cut-off time as per obturation material used or any other intended use of the device

#### **Internally Heated Tip**



- Only the working end (4 6 mm) of the tip is heated avoids need for insulating sleeves
- The non-working end is insulated and stays cool avoiding soft tissue burn



#### Class-leading Controlled Heating

- Fast heating Quick heating prevents PDL damage
- Swift cooling Swift cooling ensure GP does not follow along the heating tip

#### **Experience Ergonomics**



- 360° start / stop ring button allows holding the device in any angle comfortable for you
- Left right handed operation
- Comfortable grip, light in weight

#### **Powerful Battery**



- Highly efficient 2600mAh Lithium battery
- Direct Plug & Play In case of low battery, use the device directly with adaptor

#### **Specifications**

Heating Time	0.2 seconds for 200 °C 0.5 seconds for 250 °C
Cooling Time	2 seconds for cooling
Charging Time	2 hours to complete charge
Operating Time	4 hours of continuous operation at 250 °C








Medium (Yellow) 50/0.05 - 1 Qty

Large (Blue)

60/0.06 - 1 Qty





#### **Protocol for Use**

#### Master Cone Verification & Resizing

Select a master cone which corresponds to the master apical preparation. Check for tug-back and confirm the same with radiograph. Resize the GP by cutting it to 0.5 to 1mm short of the working length. This will prevent peri-apical extrusion of the GP

#### Selection of Fast-Pack & Plugger Tip

Select the Fast-Pack tip which reaches 4-5 mm short of the working length and give a snug fit at the tip. Mark this length using the silicon stopper. Similarly, select the hand plugger which reaches 4-5mm short of the working length

#### Gutta Percha & Sealer

Coat the master cone apically with the sealer and insert it into the canal. Brush the canal wall using strokes of GP and sealer

#### Down Pack with Fast-Pack

Remove the GP from the orifice/coronal part with the help of heated Fast-Pack tip. Plug the GP cone with pre-fit hand plugger at the orifice. Now take the heated Fast-Pack tip 4-5mm short of the working length. The bite of GP will come out of the canal along with the tip. Heated Fast-Pack tip will transfer heat to the GP and convert it into thermoplastic state. Condense the GP with the help of hand plugger to ensure the complete sealing of accessory canals and apex

#### Fill with Fast-Fill

Filling the canal using the Fast-Fill device in increments of 4mm till the canal is filled completely while condensing the warm GP with the appropriate plugger after each increment

Fast-Pack & Fast-Fill Fast-Pack & Fast-Fill 36



# UItra X

Cordless Ultrasonic Activator

#### Available Tips

Size	Length(mm)
#20 2%	18 & 21
#25 2%	18 & 21





37




Flexible Tips






## In Action



It is the process of formation of millions of gas bubbles in the irrigation solution. These gas bubbles grow in size, oscillate in different directions and explode into the liquid, producing powerful cleaning

#### Acoustic Streaming

It is a unidirectional flow of liquid produced by ultrasonic waves. Movement of tip produces special effect called acoustic turbulence which accelerates the movement in liquid inside the root canal



#### Indication

#### Smear Layer Removal

Effective in removal of smear layer and unplugging of dentinal tubules

#### Dentin Chips Removal

Ultrasonic agitation forms suspension of residual dentinal chips and tissue remnants with irrigant and facilitates their easy removal



#### Indication

#### **Biofilm Disruption**

Ultrasonic cavitation disrupts the biofilm and helps to reduce bacteria load in the canal

#### Cleaning of Complex Anatomy

Ultrasonic activation helps irrigants to flush into undetected or un-instrumented canal anatomy like isthmus, fins, lateral canals, apical ramification etc. Also helps to remove vapour lock



#### Protocol

- Prepare the root canal up to required size with the help of shaping and finishing instruments
- Choose the desired disinfectant and flushing agents for cleaning, debridement and disinfection of canals Fill the chosen irriation solution into the canal, never run Ultra X without irrigating solution
- Select the Ultra X tip of desired size and length which passively fits into canal, always keep 2 3mm short of apex Activate solution using short vertical up & down strokes for 30 seconds. Be careful to not breach apical 2mm Use the device gently and never apply pressure or wedge the tip inside the canal. If any points of resistance, then retrieve
- Replenish the irrigating solution, use suction to remove the loose debris. Repeat the cycle 4 to 5 times
- Completely dry the canal and fll it with choice of obturation material
- It is always recommended to isolate the tooth using rubber dam and use personal protective barriers like head gear, mask,



Application Based

Multiple Tips

1500 mAh

Powerful Battery





Light Weight &

Ergonomic





Cordless

Irrigation Device





Battery




4.5 Hrs of Continuous Operation



Option of Direct Plug & Play



## UltraMint & UltraMint Pro



#### **Soft Perio Power for Ultimate Comfort**


Autoclavable plastic box tips\*6, universal wrench\*1, EMS type handpiece\*1



Autoclavable plastic boxtips\*6, universal wrench\*1, Satelec type handpiece\*1



Smart touch system-ultrasensitive & waterproof





Two ways of water supply-Changing supply methods is easier than ever before



900ml stand-alone bottle with clear 20ml scale mark makes it convenient for formulating liquid.



Detachable handpiece with LEDautoclavable, ergonomic, brighter visibility





Detachable and autoclavable handpiece holder, resists cross-infection




Water control directly from handpiece






Ergonomic foot pedal











\*Satelec and EMS are trademarks of Acteon \*and EMS \*, respectively.

43 UltraMint & UltraMint Pro UltraMint Pro 44



## CuringPen

4 LEDs curing light

( Eighteeth

CuringPen

## Technical parameter:

Dimension: 24\*21\*230(mm)


Net Weight: 109g Lens size: 11mm

Power supply: >1200 times, 3 Sec Maximum light intensity



Smart controlling chip keeps light output consistent Light intensity: 1000mw/cm², 1500mw/cm², 2300mw/cm² 600 mw/cm²(detect mode)





Wavelength:385-515nm 3 Blue, 1 Violet, 10w LED



360° rotatable head



Intensity digital display




Multiple Curing Options

Ramp mode

1000 light intensity(mw/cm²)

0 5 10 15 20 Time(s)



#### Caries detection





