




**Thinner Clearer Smarter** 

### Complete imaging process within 3s Efficient diagnosis and treatment

Stable and reliable AED exposure control technology
Complete digital imaging process within 3s
Eliminating the redundant procedures in using traditional film, the images can be obtained immediately after photographing.
Get inspection results quickly






# Application in aerospace field APS CMOS technology

Efficient information access

Detailed image of low noise and high resolution





# Adopt direct-deposit Cesium Iodide (CsI) technology

Reduced signal crosstalk, more focused light output, wide adjustment range of exposure time

Low dose enables HD image acquisition, safeguarding your oral health

Ergonomically optimized Cut Corners, bringing high patient occlusal comfort

4.5mm thinnest sensor



Ergonomically optimized Cut Corners, bringing high patient occlusal comfort 4.5mm thinnest sensor allowing easier placement in the mouth








**Cut Corner** 

#### Withstand the "torture" Withstand "waves"

80000+ bending tests

80000+ bending tests
Durable and reliable
Robust CNC shell with fusion ultrasonic welding technology



#### Software

Smart sharing Rich and wonderful vision

### i-Sensor H1

**Thinner Clearer Smarter** 

