

DIAMOND BURS FG | FRESAS DIAMANTADAS FG

ROUND

ESFERICA

REF.	1010	1011	1012	1013	1014	1015	1016	1017	1018	1019	1026HL	1014HL	1016HL	1019HL	3016HL	3017HL	3018HL	1302
ISO/FIG.	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801	001/801
Ø	008	009	010	012	014	016	018	021	023	016	010	014	018	016	018	025	029	010
L1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
L2	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0	21.0	25.0	25.0	25.0	25.0	25.0	25.0	25.0	16.5
GRIT SIZE GRANULACIÓN																		

ROUND + NECK

ESFERICA CON COLLAR

REF.	1022	1023	1024	1025
ISO/FIG.	002/802	002/802	002/802	002/802
Ø	010	012	014	016
L1	2.6	2.6	2.6	3.0
L2	19.0	19.0	19.0	19.0
GRIT SIZE GRANULACIÓN				

INVERTED CONE

CONO INVERTIDO

REF.	1312	1030	1031	1032	1033	1034	1035	1036	1037	1150	1151	1153
ISO/FIG.	010/805	010/805	010/805	010/805	010/805	010/805	010/805	010/805	010/805	233/830	226/807	226/807
Ø	010	008	009	010	012	014	016	018	021	012	012	016
L1	0.9	0.9	0.9	1.0	1.2	1.4	1.6	1.8	1.8	4.0	3.5	4.0
L2	16.5	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0	21.0	21.0
GRIT SIZE GRANULACIÓN												

DOUBLE INVERTED CONE

DOBLE CONO INVERTIDO

REF.	1042	1043	1045	1046	1047
ISO/FIG.	019/806	019/806	032/813	032/813	032/813
Ø	010	012	012	014	018
L1	2.6	2.6	1.5	2.0	2.2
L2	19.0	19.0	19.0	19.0	19.0
GRIT SIZE GRANULACIÓN					

SPECIALS

ESPECIALES

REF.	2137	3138
ISO/FIG.	194/820	466/820
Ø	005	010
L1	1.0	1.0
L2	19.0	19.0
GRIT SIZE GRANULACIÓN		

BARREL

BARRIL

REF.	3180	3833
ISO/FIG.	038/811	039/811
Ø	031	034
L1	4.5	5.5
L2	19.0	21.0
GRIT SIZE GRANULACIÓN		

FLAME

LLAMA

REF.	1111	3117	3118
ISO/FIG.	243/368	243/368	243/368
Ø	014	020	023
L1	5.0	5.0	5.0
L2	19.0	19.0	19.0
GRIT SIZE GRANULACIÓN			

WHEEL

RUEDA

REF.	1051	1052	1320	3053	3054
ISO/FIG.	040/815	040/815	040/815	068/909	042/820
Ø	012	018	014	040	042
L1	0.5	0.5	0.7	1.4	0.6
L2	19.0	19.0	16.5	19.0	21.0
GRIT SIZE GRANULACIÓN					

GROSS REDUCTION

MARCADOR PROFUNDIDAD

REF.	3213	4141	4142	4217
ISO/FIG.	511/651	552/834	552/834	513/653
Ø	016	016	018	018
L1	5.0	6.0	6.0	8.5
L2	19.0	21.0	21.0	21
GRIT SIZE GRANULACIÓN				

PEAR

PERA

REF.	1147	1148	3168
ISO/FIG.	237/830	237/830	237/379
Ø	012	012	023
L1	2.5	3.5	4.0
L2	19.0	19.0	19.0
GRIT SIZE GRANULACIÓN			

FLAT END CYLINDER

CILINDRICA EXTREMIDAD PLANA


REF.	1342	1343	1090	1091	1092	1093	1094	1095	1096	1098	2173	2096	3097	3098	3099	3100	3101	4102	4103
ISO/FIG.	108/835	108/835	109/835	108/835	109/835	109/835	109/835	110/836	108/835	109/835	150/839	110/836	111/837	111/837	111/837	111/837	112/837	112/837	112/837
Ø	009	010	008	009	010	012	014	012	012	016	012	014	010	014	016	016	014	014	016
L1	2.0	3.0	4.5	3.0	4.0	4.0	4.0	6.0	3.0	4.0	0.6	6.0	8.0	8.0	8.0	8.0	10.0	10.0	10.0
L2	16.5	16.5	19.0	19.0	19.0	19.0	19.0	19.0	16.5	19.0	19.0	21.0	21.0	21.0	21.0	22.0	25.0	25.0	24.0
GRIT SIZE GRANULACIÓN																			



DIAMOND BURS FG | FRESAS DIAMANTADAS FG

CYLINDER OGIVAL END


CILINDRICA EXTREMIDAD OJIVAL



REF.	2214	2215	3215	3216	3217	3218	4220	4230
ISO/FIG.	288/ 877	288/ 877	289/ 878	289/ 878	289/ 878	289/ 878	290/ 879	290/ 879
Ø	010	012	010	012	014	016	016	018
L1	6.0	6.0	8.0	8.0	8.0	8.0	10.0	10.0
L2	21.0	21.0	22.0	22.0	22.0	22.0	24.0	24.0
GRIT SIZE GRANULACIÓN	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●

CYLINDER ROUND END


CILINDRICA EXTREMIDAD REDONDEADA



REF.	1141	2143	3145	3146	4219
ISO/FIG.	139/ 838	140/ 880	141/ 881	141/ 881	290/ 879
Ø	010	012	012	014	014
L1	4.0	6.0	8.0	8.0	10.0
L2	19.0	19.0	22.0	22.0	25.0
GRIT SIZE GRANULACIÓN	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●

TORPEDO CYLINDER


TORPEDO



REF.	2121	3122	4123	4124
ISO/FIG.	129/ 884	130/ 885	131/ 886	131/ 886
Ø	012	012	014	016
L1	6.0	8.0	10.0	10.0
L2	21.0	21.0	25.0	25.0
GRIT SIZE GRANULACIÓN	●●●●●	●●●●●	●●●●●	●●●●●

FLAT END TAPER

CONICA EXTREMIDAD PLANA



REF.	1061	1062	1063	1064	1065	1066	1332	1333	1377	2067	2068	3127	3210	3332	3069	3070	3071	4072	4200	5015	5016	
ISO/FIG.	169/ 845	170/ 845	170/ 845	170/ 845	170/ 846	170/ 846	169/ 845	169/ 845	170	172/ 847	172/ 847	170/ 845	173/ 848	172/ 847	172/ 847	173/ 848	173/ 848	173/ 848	173/ 850	173/ 878	173/ 850	
Ø	009	010	012	012	012	014	010	012	008	014	016	018	012	021	018	014	016	018	010	018	018	
L1	3.0	4.0	4.0	4.0	5.0	5.0	3.0	3.0	4.0	8.0	8.0	4.0	10.5	8.0	8.0	10.0	10.0	10.0	11.0	11.0	11.0	
L2	19.0	19.0	19.0	19.0	19.0	19.0	16.5	16.5	16.5	22.0	22.0	19.0	22.0	22.0	22.0	24.0	24.0	21.0	24.0	25.0	25.0	
GRIT SIZE GRANULACIÓN																						






DIAMOND BURS FG | FRESAS DIAMANTADAS FG

CONE ROUND END


CONICA EXTREMIDAD REDONDEADA



REF.	2131	2133	2134	2135	2136	2138	2046	3038	3131	3139	4135	4136	4137	4138	5017
ISO/FIG.	544/845	196/849	198/885	198/856	197/855	198/856	544/845	198/856	544/845	198/856	199/856	199/850	198/856	199/850	200/850
Ø	021	015	014	016	016	024	022	018	025	021	014	016	025	018	018
L1	2.3	4.0	8.0	8.0	6.0	8.0	4.0	8.0	4.0	8.0	10.0	10.0	8.0	11.0	11.5
L2	19.0	21.0	22.0	21.0	21.0	21.0	19.0	22.0	19.0	21.0	24.0	24.0	21.0	24.0	24.0
GRIT SIZE GRANULACION	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●

SAFE END


CONICA EXTREMIDAD INACTIVA



REF.	2082	3082	3083	4083	4084
ISO/FIG.	190/851	220/857	190/851	190/851	190/851
Ø	016	018	016	018	018
L1	8.0	10.0	10.0	11.0	10.0
L2	22.0	24.0	24.0	25.0	21.0
GRIT SIZE GRANULACION	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●

POINTED TAPER


TORPEDO CONICA



REF.	2112	2224	3113	3203	3205	3207	3227	3228
ISO/FIG.	297/877	297/877	298/878	299/879	299/879	299/879	298/878	298/878
Ø	012	016	012	012	014	016	016	018
L1	6.0	6.0	8.0	10.0	10.0	10.0	8.0	8.0
L2	19.0	21.0	22.0	25.0	25.0	24.0	22.0	22.0
GRIT SIZE GRANULACION	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●

NEEDLE

AGUJA



REF.	1112	1190	2200	2202	3195	3193	3200
ISO/FIG.	296/976	295/890	297/858	298/858	299/859	298/859	172/858
Ø	014	010	010	014	016	016	010
L1	5.0	3.2	7.0	8.0	10.0	8.0	8.0
L2	19.0	19.0	22.0	22.0	24.0	25.0	21.0
GRIT SIZE GRANULACION	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●

FLAME

LLAMA



REF.	4114
ISO/FIG.	249/868
Ø	016
L1	8.0
L2	22.0
GRIT SIZE GRANULACION	●●●●●

